

# Surgery Collection

Enhance your  
Product choice  
with a selection  
from Panadent

Take a look at our  
**BDIA Dental Showcase  
Special Offers** inside and  
find Panadent on **stand K03**  
Oct 9-11, 2014



Our Products, Your Success

dbg

The  
Dental  
Directory

Dentaquip

DentalSky

HENRY SCHEIN®

KENT  
EXPRESS



Trycare

RDS

Wright  
Ortodontics


To order call Panadent or your preferred dealer

t: 01689 881788 e: [info@panadent.co.uk](mailto:info@panadent.co.uk) or visit [www.panadent.co.uk](http://www.panadent.co.uk)

All Special Offers only available during the BDIA show or while stocks last. Prices subject to VAT and pp.

## CAD/CAM

### Dreve GumQuick scan




Addition-vulcanising silicone for the production of flexible and scannable gingiva masks. Optimum scannability/ readability by all commercially available dental scanners, produces excellent reproduction of detail and has a long dimensional consistency.

2 x 50ml double cartridges with 12 mixing tips, light blue & 10ml bottle

**£45.20**

### Dreve Blue Spray

Ideal for optical scans (intra-oral or of the model). Spray rich in contrast, best visibility, very long application nozzle, excellent scanning, easy to remove and washable.



50ml

**£65.70**

### Sirona Optispray

Easy to use spray for optical impressions.



50ml

**£38.00**

200ml

**£75.09**

### VITA Powder Scan Spray

A blue powder with excellent even coverage. A drying spray for admitting tooth surfaces prior to taking dental imprints.



75ml

**£89.90**

### VITA Cerec Propellant

CFC-free propellant for applying the Cerec powder to the tooth by means of the spray head.



250ml - clear gas

**£20.00**

### VITA Cerec Spray Head for Propellant

To be attached to the Cerec propellant. Application needle can be rotated by 360°



**£22.20**

### VITA Cerec powder

White contrast powder to eliminate reflection on the surface of the tooth for optical impression.

12g

**£12.30**



### VITA Cerec Liquid

Serves as a temporary adhesive for the Cerec powder. It can be completely removed with water.

6ml

**£16.80**



### VITA Ceramic Etch

Hydrofluoric acid (up to 10%) gel for etching feldspar and glass ceramic restorations. For oral use only.

6ml

**£40.30**



### VITA Simulate preparation Material

A light curing composite for the manufacture of artificial dies



(in 6 shades) for simulating the shade of the prepared tooth. It enables the dental technician to reproduce the shade more easily and accurately, since the shade of the restoration can be controlled in advance and corrected if necessary.

Preparation set

**£203.00**

## Bite registration & Impression Materials

### StoneBite®

Bite registration silicone material with an extremely high final hardness which gives an absolute precise replication of the patients' bite. One of hardest bite registration materials that breaks with clearly defined edges, easily mixes through the cannulas and has a short 80 seconds intra-oral set time.



2 X 50ml, including 12 mixing tips and 12 intra-oral tips

**£28.61**

### StoneBite® Scan

Gives perfect results with the Cerec CAD/CAM. Outstanding scanning characteristics for antagonist bite scanning. Powdering is unnecessary and it increases scan-precision with an ideal retention period inside the mouth with excellent cutting and milling properties.

2 X 50ml, including 12 mixing tips and 12 intra-oral tips

**£49.30**



# VITABLOCS® The Bedrock of Cerec for more than 25 years

## 1 Vitablocs® Mark II & Triluxe Forte

**From mill to mouth in less than 10 minutes.** Fine-particle feldspar ceramic blocks for the Cerec and inLab system. Adjust and polish intraorally, No crystallization or sintering required and less brittle than competitor porcelains offering a better finish at the incisal edges.

Mark II Available in 10 x VITA 3D shades Packs of 5	Triluxe Forte, Available in 3 x VITA 3D shades Packs of 5
Size 8, 10 & 12 Size 14	Size 12 Size 14 Size 14/14
<b>£52.40</b> <b>£60.40</b>	<b>£68.30</b> <b>£78.90</b> <b>£89.90</b>

## 2 Vita Suprinity

A NEW generation of glass ceramic material products enriched with zirconia. The result is the world's first zirconia reinforced lithium silicate ceramic. Suprinity does not need firing support, saving time and money during glazing. Suprinity ensures smoother margins and is far more easier to polish in the mouth than similar materials.

Available in 8 x VITA classical shades. Packs of 5

High Translucent	<b>£91.30</b>	Translucent	<b>£91.30</b>
<b>Sample Pack Offer: 2 Blocks with Akzent Plus LT Glaze</b>		<b>£29.99</b>	
<i>(offer only available during the BDIA show)</i>			

## 3 Vita Enamic

The first dental hybrid ceramic in the world that combines the strength of ceramic and the flex of polymer.

Available in 5 x VITA 3D shades. Packs of 5 Enamic Sets  
In High Translucent or Translucent levels

10mm, High Translucent	<b>£87.00</b>	Enamic Clinical Polishing Set	<b>£45.50</b>
14mm, Translucent	<b>£96.40</b>	Enamic Starter Set	<b>£1022.50</b>
14mm, High Translucent	<b>£101.20</b>	Enamic Staining Set	<b>£370.00</b>

**OFFER - Buy 5 x (pk of 5) Enamic or Suprinity and get 1 x Super DiGlaze diamond polishing paste and 1 x Clinical polishing set \*FREE OF CHARGE**  
**OFFER - Or Buy 5 x (pk of 5) Enamic HT and get 1 x Suprinity Trial Pack incl Akzent Glaze Spray \*FREE OF CHARGE**

## 4 Vita Akzent Plus

For simple and effective approach to stain glazing, colouring or just glazing. System comprises of stain glazing masking and glazing shades.

Available in powder, paste or spray

Powder Kit	<b>£337.00</b>	Spray Kit	<b>£427.80</b>
Paste Kit	<b>£364.20</b>		

# VITA Shade VITA Made! Expert Shade Taking with VITA

## 5 Vita Easyshade Advance 4.0

Easyshade Advanced 4.0 is the fourth generation of Easyshade technology carrying optimal shade determination from new functions for automatic measurement. With cutting edge measurement technology and durable LED components, VITA Easyshade Advance 4.0 offers precision that is twice as effective as that of the human eye.

A new "average" mode in the Easyshade Advance 4.0 allows users to measure up to 30 shades on a single tooth and arrive at an overall average shade, making the determination of a base shade on tricky anteriors even simpler. The Easyshade Advance 4.0 will automatically measure when proper placement is achieved. Shade results are displayed in VITA classical A1-D4 and VITA 3D-MASTER shades. Bleaching effects can also be displayed in accordance with the internationally recommended scale from the American Dental Association.



**Buy 5 Get 1 Free**  
(must be in the same price range)

**NEW!**



**FREE OF CHARGE**  
**Offers Available\***



**NEW!**



**Special Offer Price**  
**£995.00** (rrp £1350.00)

## Shade Guides

<b>Classical Shade guide</b>	<b>£51.20</b>
1 <b>Combined Classical &amp; bleaching attachment</b> boxed	<b>£59.40</b>
2 <b>3D Master Shade Guide</b>	<b>£51.20</b>
<b>Bleaching tabs attachment</b> for reproducing bleaching shades	<b>£9.70</b>
<b>Bleached Shade Guide</b> for monitoring bleaching	<b>£51.20</b>
3 <b>Linear Shade Guide</b>	<b>£98.90</b>
The ease of Classical shade guide with the accuracy of VITA 3D Master System. The ultimate tool for finding the right shade, every time.	
<b>Special Offer Price</b>	<b>£89.00</b>
<b>Vitablocs Guide 3D Master</b>	<b>£51.20</b>
Shade tabs are made of original Vitabloc Mark II ceramic for a quick and accurate Vitablocs shade selection	
<b>Individual Shade Tabs</b> are all available for replacement,	from <b>£3.50ea</b>

## **NEW!** RiteLite 2 LED

RiteLite LED shade matching light, affordable, cordless and maintenance free. RiteLite provides a neutral light environment around the restorative area.

With polarizing filter - 12 LED's includes one neutral colour pad with 25 sheets	<b>£295.00</b>
Without polarizing filter - 12 LED's includes one neutral colour pad with 25 sheets	<b>£221.00</b>

## VITA Teeth

### 4 Vita Physiodens (High end)

Precise shade selection, perfect setup, with broad occlusal surface reliefs, perfect closure of interdental spaces and high quality material guarantee perfect results in every case. You save time, win customers and prestige. These teeth are unique products of high-tech craftsmanship, each manually layered in a multi-layer technique.

Anteriors set of 6	<b>£21.20</b>	Posteriors set of 8	<b>£19.30</b>
--------------------	---------------	---------------------	---------------

### 5 Vita Vitapan (Mid to High end)

A timeless classic of acrylic teeth made of synthetic resin material of precisely matched contents of inorganic fillers with gradient polymers.

Anteriors only, set of 6	<b>£17.60</b>	Posteriors set of 8	<b>£15.90</b>
--------------------------	---------------	---------------------	---------------

### 6 Vita MFT (Independent)


Multi Functional Teeth in harmony with Nature MFT gives a play of colours, rich in nuances, that catch light reflection, comfortable wearing which supports facial harmony, easy setting with Easy Centric where the surfaces are designed in such a way, almost automatically, to achieve the perfect centric.

Anteriors set of 6	<b>£8.95</b>	Posteriors set of 8	<b>£8.95</b>
--------------------	--------------	---------------------	--------------

### 7 Vita Lingoform (Universal posterior mould)

High quality denture tooth for all situations with a lingualised occlusion concept to reduce lateral space loading. Can also minimise destructive lateral forces on implant support dentures.

Posteriors set of 8	<b>£16.50</b>
---------------------	---------------



### Special Offer:

Choose 24 Anterior sets in either range and get a metal display frame **FREE OF CHARGE**

## ORAL SCREENING, HAIRLINE CRACKS/CARIES DETECTION

### Microlux/DL & Microlux/Transilluminator

"Addent has given us a handy little tool, no bigger than a dental handpiece, with a myriad of tips and supplies that assist us in diagnosing everything from initial caries lesions, micro-fractures, through to some abnormal changes in soft tissues" **David Winkler DDS**, Castle View Dental, Windsor

### Oral Cancer Screening

Could not be easier with the aid of a second pair of eyes, the DL. Affordable enough to use for every adult patient's clinical oral examination. Early identification of at risk patient's is key to reducing incidence of oral cancer. Use the DL light technique as your HIGHLIGHTER to help see abnormalities clearly and therefore earlier.



#### Offer:

**Microlux handpiece with choice of DL (Oral Cancer Screening attachment) or Transilluminator or proximal caries attachment:**

**£195.00**  
Save **£24.00**

### Other Available Attachments

#### Transilluminator

A high quality image grade glass light guide designed to focus the LED light source at the output. For the detection of hairline cracks, caries and calculus for general dentistry, prosthodontics, endodontics, periodontics, dental hygienists, oral surgery, orthodontics and pedodontics. Also ideal for pregnant patients and children.

**Transilluminator 2mm £98.00**  
micro black glass light guide attachment

**Transilluminator 3mm £98.00**  
standard black glass light guide attachment

**Metal attachment £98.00**  
for the Perio probe (w/5, 7 & 9mm light markers), Endo light (1mm) or Proximal caries (0.75mm) light fibres

**Refill box £25.00**  
25 x dbl sided fibres for either the perio probe, endo light or proximal caries

**Lighted Mirror attachment £67.00**  
box of 5

**Microlux Handpiece £125.00**  
housing a no-heat LED light source with 3 x N 1.5v batteries

All attachments are autoclavable. Light fibres are single use.

**Protective sleeves £25.00**  
for DL & Transilluminator box of 250

**Pro-C (30ml oral cleanser) £2.69**

**Little blue presentation box £15.80**



**Full kit offer £630.00**  
Save **£70.00**

*"The early detection of mouth cancer is the most important service dentists can offer their patients. Microlux/dl helps with this. It's simple mode of use can be easily learned by all clinical team members allowing the formulation of a cohesive screening protocol across the practice."* Mr Philip Lewis BDS, Avenue Road Dental Practice, Isle of Wight, an ambassador for the Mouth Cancer Foundation



### ORABLU - Toluidine Blue



Oral lesion marking system:

**For Dentists** - Use ORABLU to refer a patient with a photo of the stained lesion/patch.

**For Specialists** - use ORABLU to mark the lesion for biopsy. A very helpful stain when the biopsy position is not absolutely clear.

Meets the requirements of ARTICLE 14 OF THE Medical Devices Directive, 93/42/EEC Is designed to be used after a lesion/patch identification with suspected pathology. ORABLU is not diagnostic. A biopsy alone is diagnostic.

Specific to oral cancer see ORABU articles at [www.addent.com/articles](http://www.addent.com/articles).

ORABLU™ Oral Lesion Marking System consists of one bottle net 3.0ml containing 0.5% Tolonium Chloride, also known as Toluidine Blue, one bottle net 6.0ml containing 1% acetic acid solution (cleanser), and 30 individually packaged applicator swabs plus full instructions for use.

1 box for up to 10 applications **£72.50**

### Need help?



Call our friendly customer service team on **01689 881788** or email us [info@panadent.co.uk](mailto:info@panadent.co.uk) [www.panadent.co.uk](http://www.panadent.co.uk)

## Bite registration & Impression Materials

### Regofix Transparent

Transparent bite registration material, practically free of filler material and with high final hardness and a quick recovery from deformation. Easy to cut not opaque to x-rays.

2 X 50ml double cartridges  
incl 12 mixing tips, green



**£25.50**

### Fresh® Putty

A popular addition - vulcanising impression silicone. Non-sticky and easily kneadable. Suitable for sandwich; correction, double mixing technique and valve border impressions. Mint taste

2 x 250ml incl. dosing spoons  
2 x 450ml incl. dosing spoons



**£30.91**  
**£47.40**

### Fresh® Mono

A very precise addition-vulcanising impression silicone. Perfect impression taking with only one material.

2 x 50ml, incl. 6 mixing tips, green



**£30.90**

### Fresh® Clear

Clear, transparent, addition-vulcanising, monophasic impression silicone. Ensures a perfect control already during the impression-taking. Universally applicable and particularly suitable for the implant impression-taking.

2 x 50ml incl. 6 mixing tips, green



**£45.50**

### Zerosil Soft

Easy and fast addition... impressions. Disinfectable and compatible with plasters. Can be re-poured several times.

2 x 500ml tubes



**£46.60**

## Mouthguard & Tray Production

### Dreve Druformat Scan

High quality thermo-forming unit with the latest material recognition software. Automatic scanning device for blanks recognition. High operating pressure of 6 bar for all professional applications. Design colours: Red or Silver



**Unbeatable Offer**  
**£2100.00**

**With £400 worth of FREE goods (including the Funky Tool)**

### Dreve Drufosmart Scan

Pressure forming unit with automatic material recognition. Variably adjustable operating pressure of 2-4 bars. Suitable for blanks between 0.1mm and 5mm thickness. Design colours: Red or Silver



**Unbeatable Offer**  
**£1420.00**

### Dreve Funky Tool

An optional Tool compatible with the Druformat Scan. To recycle and create customised blanks. Mix and match colours for totally unique patterns.



**Unbeatable Offer**  
**£299.00**

### Vacformat U

Robust and compact (155x365x180mm-6.1kg) vacuum thermoforming unit suitable for all blanks with a thickness of 0.1 to 4mm.

Integrated vacuum pump, therefore only a connection to mains supply is required. Vacuum power: 0,75 bar. Steel granules are recommended to be used with this unit.



**Unbeatable Offer**  
**£815.30**

## Selection of Dreve High Quality Blanks

For positioners, retainers, bruxism splints and splint therapy.

### Drufosoft Transparent

Made of EVA. 82 shore A. High recovery from deformation.

1.5mm	20 pieces	<b>£23.00</b>
2mm	15 pieces	<b>£15.70</b>
3mm	10 pieces	<b>£23.50</b>
4mm	15 pieces	<b>£36.50</b>
5mm	15 pieces	<b>£29.00</b>



### Drufosoft Colour



3mm	Any colour, pack of 3	<b>£14.40</b>
4mm	Any colour, pack of 10	<b>£32.50</b>

### Drufosoft Bicolour

3mm	3 pieces black/yellow/black white/blue/white red/blue/red	<b>£15.50</b>
-----	--	---------------



### Drufosoft Tricolour

3mm	3 pieces black/red/yellow blue/white/red red/green/blue	<b>£18.00</b>
-----	--	---------------



### Drufosoft Pro

3mm	10 pieces transparent, made of EVA Hard-flexible, 95 shore A. High recovery from deformation	<b>£17.50</b>
-----	--	---------------



### Biobleach

For bleaching trays.

<b>Biobleach Hard</b>	Polyethylene 0.5mm 25 pieces	<b>£22.00</b>
-----------------------	---------------------------------	---------------

<b>Biobleach Soft</b>	Ethylvinyl Acetate 1mm 20 pieces	<b>£13.50</b>
-----------------------	-------------------------------------	---------------



### Biolon

Made of polyethylene. Suitable for both Vacuum and pressure forming. Clear transparent, for hard occlusion splints. Base for occlusion splints, base for orthodontic devices.

Positioners/ retainers (0.75mm). X-ray splints. Temporaries.

0.5mm	25 pieces	<b>£15.50</b>
0.75mm	20 pieces	<b>£16.50</b>
1mm	25 pieces	<b>£23.50</b>



### Kombiplast

For occlusion splints and more comfortable wearing. Hard surface and soft inner side. Transparent.

2mm	10 pieces	<b>£22.50</b>
3mm	10 pieces	<b>£29.50</b>
4mm	10 pieces	<b>£31.50</b>

## Dreve Must Have Products

For positioners, retainers, bruxism splints and splint therapy.

### MaxiFresh

25 cleaning and care tablets for all dental splints and sports mouth guards. FOR HOME USE – Abrasion free and highly recommended. Daily cleaning prolongs longevity & ensure better hygiene & function.



Box of 25 with cleaning bowl

**£9.51**

Box of 24 refill tablets

**£7.01**

### Dreve

#### EyeVolution

Perfect for light curing Enamic (3 minutes). Highly efficient and a most attractive design housing LED light waves for ultimate cold-light polymerization.



**Unbeatable Offer**  
**£489.00**

### DocSnoreNix<sup>®</sup> Anti-snoring Device

**VERY GOOD VALUE**

DocSnoreNix Starter Set  
Containing:  
Pair of spring pliers, 1 x Set of Protrusion Springs for Patient, 2 x Kombiplast 3mm thermoforming blank Protrusion Calibration, Lightdon Gel 20ml Transparent, Lightdon Bond 10ml bottle UltraTrim Medium, Accessories, Manual.



**Introductory Kit**  
**£129.50**


**Refill Pack**  
**£80.50**

## Thomas Endodontic Instruments



### Thomas T-Flex Files

Unique Stopper System safeguards puncturing the apex. Ideal for deeply curved RCT treatments and procedures where precision and reliability is imperative. These files are fitted with the adjustable stopper system delivering complete flexibility when measuring length to the apex and safeguarding against apex puncture.



#### Hand Reamers

Assortments - Length (mm) 21, 25 and 28:  
 15/40: **£8.87**    45/80: **£10.21**    90/140: **£12.20**  
 Individual pack:  
 6-10: **£8.73**    15-50: **£9.08**    55-80: **£10.38**  
 90/140: **£15.40**

#### K Hand files

Assortments - Length (mm) 21, 25 and 28:  
 15/40: **£9.25**    45/80: **£10.38**    90/140: **£12.97**  
 Individual pack:  
 6-10: **£8.90**    15-50: **£9.25**    55-80: **£10.38**  
 90/140: **£12.97**

#### H Hand files

Assortments - Length (mm) 21, 25 and 28:  
 15/40: **£9.09**    45/80: **£12.02**    90/140: **£12.02**  
 Individual pack:  
 8-10: **£9.04**    15-50: **£9.09**    55-140: **£12.02**

#### Superflex H files

Assortments - Length (mm) 21, 25:  
 15/40: **£9.89**    45/80: **£11.66**    90/140: **£12.02**  
 Individual pack:  
 10: **£9.85**    15-50: **£9.89**    55-80: **£11.66**

#### Engine RA Reamers

Assortments - Length (mm) 21, 25 and 28:  
 15/40: **£10.77**    45/80: **£12.38**    90/140: **£14.24**  
 Individual pack:  
 10: **£10.37**    15-50: **£10.77**    55-80: **£12.38**

#### Dynamatic RA Files

Assortments - Length (mm) 21, 25:  
 15/40: **£9.71**    45/80: **£11.11**    90/140: **£14.24**  
 Individual pack:  
 10: **£8.63**    15-50: **£9.71**    55-80: **£11.11**

**Special offer:**  
 Buy 5 packs in the same price range  
 Get 6th Pack FREE OF CHARGE

### Standard Paste Filler

Assortments - Length (mm) 21, 25 and 29: 25/40: **£12.37**  
 Individual pack: 025, 030, 035, 040: **£12.37**


### Special Paste Filler

Assortments - Length (mm) 21, 25: 25/40: **£13.73**  
 Individual pack: 025, 030, 035, 040: **£13.73**

## Universal Post Remover (U.P.R)

designed by Pierre Machtou

A best seller for many years. Removing posts from root canals prior to endodontic treatment is often a challenge, but the Thomas U.P.R makes light work of this task and keeps the procedure simple and expedient. The kit comes complete with all instrumentation required for nearly all cases.



Stainless box and instruments for cleaning and sterilization without rust.

Box comprising of:  
 1FG Diamond Bur; 1 Pointeau drill;  
 1 Peeso Reamer n2; 2 sets of 4 trephines (different diameters); 2 sets of 4 Taps counter clock wise; 1 set of metallic and silicone washers (different diameters); 1 Plier

**Special Offer Price**  
**£700.00** RRP **£810.32**

## NEW! Endopener

The NEW high cutting power, rotary nickel-titanium instrument, is designed to enlarge the root canal orifices without a screwing effect. With low stress on the instrument and the root canal, through a continuous rotation and reciprocating motion, it is extremely efficient.




Pack of 4

**£28.00**

## NEW! NiTiShapers

NiTi Canal Tapered Shapers (Nickel Titanium) provide a logical approach to maximise speedy & effective canal opening with clearly marked length markers & a high degree of flexibility to get around even the most curved root canals.



#### Start-Up Bundles:

1. Pack of 4 Endopeners rrp **£28.00**
2. Assortment pack of NiTiShapers T1, T2, T3 - rrp **£38.00**
3. Assortment of K-Hand Files with Stopper Handle, sizes 14-40, Length 25mm - rrp **£10.38**

**Special Start-Up Offer**  
**£55.00**

**NiTiShaper Intro Offer**  
**£32.500**





## inSafe Needle disposal system

inSafe:

### For All Dental Surgeries

Meets directive 2010/32/Eu legal requirements for all dental surgeries to have a safe disposal of sharp procedure in place.

Introductory pack contains: inSafe (not InSafe) syringe, inSafe sharp container with base and 250 inSafe needle adaptors.

Introductory Pack contains: InSafe Syringe, 2.2 ml, O Ring Aspirating + inSafe Sharps Container + inSafe Sharps Container Base + inSafe Needle Adaptor, Imperial (Yellow) 250 pcs.



**Introductory Special Offer**  
**£67.00**

**Additional syringes as part of the same purchase**  
**£45.00 only**

## Peri Implantitis

### Tigran No2

### Peribrush



For treating ZPeri-Implantitis **£57.00**

## Prime-Dent Desensitiser

Contains hema and benzalkonium chloride. Fast, effective way of sealing off dentinal tubules by combining with proteins and forming a physiologic barrier of protection. Benzalkonium chloride adds an anti-microbial action to eliminate bacterial. Used as a normal part of resin bonding. Also effective for cervical sensitivity, root planning and scaling.

7ml

**£24.99**



## Compoglaze

Has a lustrous effect on all resin materials (composite, compomers, acrylics, etc). Dedicated diamond grit size specifically for resins with rapid chairside results.

**£29.71**



## Super Diglaze

To re-glaze porcelain intra orally in just 30 seconds. An ideal polishing paste for bringing back lab gloss finish crowns, bridges, etc that have been adjusted in the mouth. It adheres beautifully to rotary mop or felt

**£29.71 ea**



## Gutta Percha Paper Points

Uniformly hand rolled using the finest grade raw materials and uncalenared paper stock for optimum performance. All are ISO colour-coded for ease of selection in a wide range of sizes.

Box of 120

**£9.65**

(please state required size)



## High Quality Implant System

otmedical®

### OT F1-F3 Implants

Made in Germany

OT Medical offer a complete implantology concept which includes high quality and innovative products.

#### New 4plus6Line for OT-F2 implant system

Based on the "Multi-Unit Treatment Concept", the users of the self-tapping OT-F2 screw implant benefit from the new "4plus6Line" abutments. Which offer the possibility to restore an edentulous mandible or maxilla with a fixed denture on only 4 to 6 implants in just one visit (immediate loading).

The posterior implants inserted in 30° angulation allow the use of greater length implants (preferably 16 mm), thereby making optimal use of the available bone. Augmentation to increase the vertical alveolar height, a sinus lift or nerve lateralization can be completely avoided.

The 4plus6Line abutments are available in 0° version as well as in angulations of 17° and 30°, each in socket heights (GH – gingival height) of 1.50 and 3.00 mm. Supplied in a sterile packaging, they can be immediately fixed to the inserted implants during surgery. The abutments contain the mounted color-coded final screw.

See the full OT-F1, F2 F3 and Bicortical range on Panadent website. [www.panadent.co.uk](http://www.panadent.co.uk)



**OT-F1**  
3.80 x 10mm  
Nanopost  
**£165.00**



**OT-F2**  
3.40 x 8mm  
Nanopost  
**£165.00**



**OT-F3**  
5 x 5mm OT3  
Press-fit Implant  
**£198.00**

## Panadent Essentials

### Minamint Original

The original highest quality Super Concentrated chairside mouthrinse. Made in England for over 20 years, now with added fluoride.

90ml **super concentrated**  
**£17.72**

Beware of imitations!



### **NEW!** Minamint Peroxide Whitening Mouthwash

Ideal for patients wishing to bring the Minamint experience home. For brightening their smile, freshening their mouth, maintaining tooth whitening following whitening procedure, implant cleanliness and general oral hygiene.

Minamint Pre-brush 2% Peroxide Whitening mouthwash. One bottle will last for one month if patient rinses twice daily.



#### Offer 1

Buy 12 x 500ml minamint whitening mouthwash at **£4.00** each  
Get 1 FREE Minamint Original

#### Offer 2

Buy 24 x 500ml minamint whitening mouthwash at **£3.90** each  
Get 3 FREE Minamint Original

#### Offer 3

Buy 48 x 500ml minamint whitening mouthwash at **£3.70** each  
Get 6 FREE Minamint Original

#### Offer 4

Buy 60 x 500ml minamint whitening mouthwash at **£3.50** each  
Get 12 FREE Minamint Original

Buy 2 Minamint Original or 2 Minamint Zinc and Get 1 Minamint Peroxide Whitening Mouthwash FREE

## A touch of Gold - Midas 'Non stick composite instruments set'

Made in Germany by Sword

Buy any 5 of the following instruments for only **£155** and receive a **Sword Frame** worth **£30.00** Free of charge

includes:

- Heideman Flexible Spatula (4086/0)
- Filling instrument for Composite (4086/101)
- Filling instrument for Composite (4086/107)
- Carving instrument for Composite (4086/122)
- Modelling Instrument for Composite (4086/123)



4086/0    4086/101    4086/107    4086/122    4086/123

## Calset Composite Heater



The luxury of Composite warming

- Base heater with 230 volts AC power supply **£257.70**
- Calset Heater, 230V AC power supply with 7 hole syringed Composite Tray **£382.00**
- Calset Heater, 230V AC power supply with Tri-Tray Composite tray **£412.00**
- For CoMax, Caulk, Kerr (Original Model), SDI and Compo Jet also holds 2 finishing instruments. Another top (not shown, is supplied for Centrix, 3M, Venus, Kerr (newer model) Guns.
- Calset Heater, 230V AC Power Supply with standard Composite Compule Tray **£374.00**

Mix & Match Trays for the Calset heater base

- Anesthetic Tray (holds 3 carpules) **£150.00**
- 7 hole Syringed Composite Tray **£150.00**
- Standard Composite Compule Tray **£150.00**
- Tri-Tray **£190.00**
- Dispenser Tray **£190.00**
- Calset Porcelain Veneer Tray **£150.00**

# Disinfection and Decontamination

## OROCLEAN (OCC) Switzerland

Highest quality Swiss cleaning & disinfection products with unprecedented material compatibility.

Ask for your **FREE Dentiro Zero Sample...**

and free OCC Infection Control hardcover manual/product catalogue with valuable tips and tricks for effective infection control, when placing an order.



**FREE Beach Bag**  
with all OCC orders over £90

## Isorapid Spray

Isorapid Spray for medical devices/ surfaces (equipment, installations & instruments), Aldehyde Free, extremely low alcohol content of less than 50%, fast acting, broad spectrum disinfectant effective against noroviruses

1 litre bottle **£9.76**  
sprayer **£1.35**  
5 litre **£45.76**



Can also be used to increase the saturation of dentiro wipes or make your own pre-moistened wipes with Isorapid dry XXL wipes.

## Isorapid XXL Dry Wipes

Pack of 70 **£7.95**

## Isorapid Floor

1 litre **£18.75**

**Buy 5 litres at £46.50**  
Get 1 extra litre and  
120 XXL Dry Wipes FREE



## Oromed Gel

hydro-alcoholic with a unique skin protection system for surgical and hygienic hand & forearm disinfection.

1 litre bottle **£19.38**



## Dentiro Classic

120 wipes for surfaces and medical devices  
14.5 x 20cm **£5.45**



## Orocid-Multisept Plus

For medical instruments and burs

1 litre bottle **£18.02**  
Dispenser **£6.06**

## Aseptoprint

For dental impression material

1 litre bottle concentrate **£50.23**  
50 litres dispenser **£6.06**  
1litre spray **£13.15**



## Dentiro Zero

New Alcohol free, for surfaces & non-invasive medical devices (equipment, installations & instruments)

1 litre **£10.50**  
5 litre **£46.50**



## Rusnon Medical

Anti-corrosive and cleaning disinfectant formulated for endodontic instruments and diamond burs. Use neat or diluted. 50% concentrated. Contains 20% chlorhexide. Triple action solution for cleaning, protecting and sanitising rotary and endodontic instruments.

1200ml bottle **£26.83**



## Omnisan Plus

For disinfection of aspirators and spittoons. Kill all microorganisms in 4 minutes. Non foaming and anti corrosive solution which contains sodium dichloroloscyanurate

60 Effervescent Tablets  
for 60 litres. **£12.93**



## **NEW!** Orolin Intra

Efficient and Cost-effective Cleaning and Disinfection of fast and slow Hand Pieces (popular brands).

### HTM 01-05 Compliant

Estimated cost of disinfecting and oiling each hand-piece -10-12p

**Introductory Offer**  
**£48.70** for slow handpiece



## Adaptors for Fast hand piece

Kavo **£15.50**  
NSK **£15.50**  
Bein Air **£22.99**  
Sirona **£19.51**  
W&H **£19.00**  
Siemens **£24.00**



**NEW!**

## Sirona Apollo DI - Intra Oral Scanner

Apollo DI has a small, easy access video camera head, offering precise, high quality imagery quickly and accurately. The images are seamlessly forwarded in open STL format via Sirona's internet portal, 'Sirona Connect', to the technician to process. There are no hidden licence charges and software updates are free of charge. Milling is performed in the Laboratory.

**Apollo DI**  
**£13,650.00**




**NEW!**

## Bioart Articulator 4000S Semi-Adjustable Articulator

The most traditional model provides the key features of a semi-adjustable articulator, including intercondylar distance adjustment via micrometric expanding spindle. The Elite facebow brings an articulated and fixation system for the fork that allows a faster registration process and more comfort for the Patient.



**Bioart Articulator 4000S**  
**£275.00** with NEW Elite Face bow  
and jig transfer in carry case



Our Products, Your Success

To order call Panadent or your preferred dealer

t: 01689 881788 e: [info@panadent.co.uk](mailto:info@panadent.co.uk) or visit [www.panadent.co.uk](http://www.panadent.co.uk)

All Special Offers only available during the BDIA show or while stocks last. Prices subject to VAT and pp.